



### Staining Station Instructions:

- Make arrangements for training with Shreyas prior to using the staining station for the first time.
- Please ensure that during your orientation, you understand both the staining procedure and how to maintain the staining set-up.
- If back up/replacement of solution/s is needed, it is best to do that just before you start staining. **Talk to Shreyas about this.**
- **You must book a slot for staining using the booking calendar at least 24 H in advance.**
- Make sure that the fume adsorber is turned on before you start the process of staining.
- You can stain up to 24 slides at a time using the slide holder.
- Please record the number of slides and racks of slides (1-24 slides in 1 rack) that have been stained in the logbook found above the weighing scales. This information is needed to keep track of usage and check status of reagents.
- To prevent excessive carryover of reagent from one staining jar to the next, please drain back as much reagent as possible.
- For staining fees per slide, refer to cost recovery pricelist on the histology core facility website.
- **The fume adsorber must be left turned on for at least 3 hours after staining. If working after hours or on weekends, leave the fume adsorber turned on overnight or over the weekend.**

### Coverslipping Station Instructions:

- The cost of coverslipping is included in the overall cost of staining for one slide.
- First time users must seek training for coverslipping from Shreyas.
- Make sure that the fume adsorber is turned on before you start the process of coverslipping.

- If you are running low on mounting media, ask Shreyas or alternatively transfer the Acrytol mounting media (Leica Inc.). This can be found in the cabinet above the coverslipping station with a label on the door indicating mounting media.
- We provide a variety of sizes of coverslips. **Please use the long coverslips (e.g., 24x50 mm or 24x60 mm) only if you have those many sections on a slide requiring coverslipping.**
- **AFTER COVERSLIPPING, SLIDES MUST BE PLACED IN THE FUMEHOOD OVERNIGHT FOR DRYING OUT! THE MOUNTING MEDIA CONTAINS XYLENE AND CAN RELEASE NOXIOUS FUMES - NO EXCEPTIONS!!!!**
- If you run out of coverslips, you can grab a box of coverslips from the drawer labelled coverslip (below IHC-2 work bench and opposite to the microwave).

## Alcian Blue and Nuclear Fast Red Staining Protocol for Paraffin Embedded Tissues:

Slides are warmed to room temperature for at least 1 hour to get them to room temperature. After slides are thawed at room temperature, they are baked at 60 °C for at least one hour or overnight prior to the next step i.e., deparaffinization. The slide baking can be conducted in Oven#1 (close to the main lab entrance).

1. 3 Minutes Xylene – 1
2. 3 Minutes Xylene – 2
3. 1 Minute Xylene / Absolute Alcohol
4. 1 Minute Absolute Alcohol – 1
5. 1 Minute Absolute Alcohol – 2
6. 1 Minute 95% Ethanol
7. 1 Minute Tap water – WHITE BUCKET
8. Rinse Distilled Water – BLUE BUCKET
9. 30 Minutes Alcian Blue solution
10. Wash Tap Water (2-3 changes) – WHITE BUCKET
11. 2 Minutes Running Tap Water – WHITE BUCKET
12. Rinse Distilled Water – BLUE BUCKET
13. 5 Minutes Nuclear Fast Red Solution
14. 1 Minute Running Tap Water – WHITE BUCKET
15. 3 Minutes 95% Ethanol
16. 3 Minutes Absolute Alcohol – 1
17. 3 Minutes Absolute Alcohol – 2
18. 1 Minute Absolute Alcohol / Xylene
19. 1 Minute Xylene – 1
20. 1 Minute Xylene – 2
21. 1 Minute Xylene – 3

Leave slides in xylene until ready to coverslip; don't let them dry out!

Coverslip slides in a permanent mounting medium in the Coverslipping station only! **AFTER COVERSLIPPING, SLIDES MUST BE PLACED IN THE FUMEHOOD OVERNIGHT FOR DRYING OUT! THE MOUNTING MEDIA CONTAINS XYLENE AND CAN RELEASE NOXIOUS FUMES - NO EXCEPTIONS!!!!**

**Results:**

Strongly acidic sulfated mucosubstances (goblet cells of intestine) ----- blue

Nuclei ----- pink to red

Cytoplasm ----- pale pink

**Reagents:**

3% Acetic Acid Solution:

Glacial acetic acid ----- 3 ml

Distilled water ----- 97 ml

Alcian Blue Solution (pH 2.5):

Alcian blue, 8GX ----- 1 g

Acetic acid, 3% solution ----- 100 ml

Mix well and adjust pH to 2.5 using acetic acid.

0.1% Nuclear Fast Red Solution:

Nuclear fast red ----- 0.1 g

Aluminum sulfate ----- 5 g

Distilled water ----- 100 ml

Dissolve aluminum sulfate in water. Add nuclear fast red and slowly heat to boil and cool. Filter and add a grain of thymol as a preservative.

**Revisions:**

October 2023- The protocol was updated from the references below for the histology core facility by Shreyas Jois and Adi Manek with facility specifics.

October 2023

## References:

1. <https://healthsciences.usask.ca/facility-services/Histology/alcian-blue-and-alizarin-red-s-ffpe-tissues-version-0.1-july-2022.pdf>
2. [http://www.ihcworld.com/\\_protocols/special\\_stains/alcian\\_blue.htm](http://www.ihcworld.com/_protocols/special_stains/alcian_blue.htm)
3. [https://www.ihcworld.com/\\_protocols/counterstain\\_solutions/nuclear\\_fast\\_red.htm](https://www.ihcworld.com/_protocols/counterstain_solutions/nuclear_fast_red.htm)